

Quick Notes on Dotter usage

Files used with Dotter

DNA1_sequence_for_dotter
DNA2_sequence_for_dotter
PEPTIDE1_sequence_for_dotter
PEPTIDE2_sequence_for_dotter

Commands used

% dotter DNA1_sequence_for_dotter DNA2_sequence_for_dotter
% dotter PEPTIDE1_sequence_for_dotter PEPTIDE2_sequence_for_dotter
% dotter DNA1_sequence_for_dotter PEPTIDE1_sequence_for_dotter

Results

The results are graphical, but if you need to save results, just issue the command:

dotter -b <results_file_name> <horizontal_sequence> <vertical_sequence>